IMMUNOBIO

4 in 1 SARS-CoV-2 Antigen Rapid Test

BfArM Test-ID: AT228/20 (naso-/oropharyngeal) | AT150/21 (anterior nasal) | AT149/21 (Speichel)

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Manufacturer:

Name: Hangzhou Immuno Biotech Co.,Ltd.

Address: Room 401, 4F, Building 2, No.28, No.3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China.

Tel: 0571-85368996

Email: info@jgbiotech.com

We, the manufacturer, here with declare that the product(s)

Product Name	SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag)	Specification	40 T/kit, 20 T/kit, 10 T/kit, 1 T/kit.
Intended Use	The SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) is a rapid chromatographic immunoassay for the qualitative detection of novel coronavirus SARS-CoV-2 in human throat secretions, nasal secretions and saliva specimen.		
Classification	Others		

Conformity Assessment Route : IVDD98/79/EC Annex III.

Applicable Standards:

ISO 13485:2016

ISO 14971:2019

EN ISO 18113-1:2011

EN ISO 18113-2:2011

EN ISO 18113-3:2011

EN 13641:2002

ISO 15223-1:2016

EN 13612:2002

ISO 23640:2015

EN 62366-1:2015

We, the manufacturer, here declare with sole responsibility that our product/s mentioned above meet/s the provisions of the Directive 98/79/EC of the European Parliament and of the Council on In-Vitro Diagnostic Medical Devices.

We agree to develop, implement and maintain a documented post-production monitoring process.

Name of General Manager	Xiaoguang Li
Signature	
Date	Nov 30, 2020
Place	Hangzhou, China.
Seal (Manufacturer)	

This is to certify that the Quality Management System of

Hangzhou Immuno Biotech Co., Ltd.

Unified Social Credit Code: 91330101092042370C

Operation Address: Room 401, 4F, Building 2, No.28, No.3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China

Registered Address: Room 401, 4F, Building 2, No.28, No.3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China

applicable to

R&D, Production and Sales of Coronavirus COVID-19 IgG/IgM Antibody Rapid Test, SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag), SARS-CoV-2 Neutralizing Antibody Rapid Test (COVID-19 Ab)(Export to EU)

has been assessed and registered by NQA against the provisions of

ISO 13485: 2016

This registration is subject to the company maintaining a quality management system, to the above standard, which will be monitored by NQA.

Certified Clients shall accept regular surveillance assessments, the validity of certificates shall be maintained for the positive result of audit.

The information of this certificate can be checked on CNCA's website (www.cnca.gov.cn)

SNQA's website: www.snqa.com.cn

Managing Director

Certificate Number **47938**

Date: 27 July 2020
Reissue Date: 13 January 2021
Valid Until: 27 July 2023

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Hangzhou Immuno Biotech Co.,Ltd.

Hinzufügen: Room 401, 4F, Building 2, No. 28, No. 3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China

Tel: 0571-85368996

E-mail: info@jgbiotech.com


IMMUNOBIO

SARS-COV-2 Antigen Rapid Test COVID-19 Ag

SARS-CoV-2 Antigen Rapid Test (COVID-19) is a rapid chromatographic immunoassay for the qualitative detection of novel coronavirus SARS-CoV-2 in human nasal cavity.

Hangzhou Immuno Biotech Co.,Ltd.
Room401, 4F, Building 2, No.28, No.3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China

CONTENT OF THE KIT

PRECAUTIONS FOR USE


- For in vitro diagnostic use only. Do not use after expiration date.
- The Test strip should remain in the sealed pouch until use.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- The used Test Strip should be discarded according to national, state and local regulations.
- Humidity and temperature can adversely affect results.

INTERPRETATION OF RESULTS



- Positive (+): Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the T line region.
*NOTE: The intensity of the colour in the test line regions may vary depending on the concentration of SARS-CoV-2 present in the specimen. Therefore, any shade of colour in the test line region should be considered positive and recorded as such.
- Negative (-): One coloured line appears in the control line region (C). No line appears in the T line region.
- Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test.
If the problem persists, discontinue using the test kit immediately and contact your local distributor.

- PLEASE READ THE INSTRUCTIONS (and the information in the additional leaflet) THOROUGHLY BEFORE STARTING THE TEST PROCEDURE
- Have a watch or a timer (only equipment not provided).
- Take all the elements out of the box, place them on a flat and clean surface and identify each element of the kit (see content on the left).

SPECIMEN COLLECTION

TEST PROCEDURE


IMMUNOBIO

SARS-COV-2 Antigen Rapid Test COVID-19 Ag

SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) is a rapid test for the qualitative detection of antigens to the novel coronavirus SARS-CoV-2 in sputum.

Hangzhou Immuno Biotech Co., Ltd.
Room 401, 4F, Building 2, No. 28, No. 3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China

CONTENT OF THE KIT

PRECAUTIONS FOR USE


- For in vitro diagnostic use only. Do not use after expiration date.
- The Test strip should remain in the sealed pouch until use.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- The used Test Strip should be discarded according to national, state and local regulations.
- Humidity and temperature can adversely affect results.

INTERPRETATION OF RESULTS

- Positive (+): Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the T line region.
- *NOTE: The intensity of the colour in the test line regions may vary depending on the concentration of SARS-CoV-2 present in the specimen. Therefore, any shade of colour in the test line region should be considered positive and recorded as such.
- Negative (-): One coloured line appears in the control line region (C). No line appears in the T line region.
- Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

- PLEASE READ THE INSTRUCTIONS (and the information in the additional leaflet) THOROUGHLY BEFORE STARTING THE TEST PROCEDURE
- Have a watch or a timer (only equipment not provided).
- Take all the elements out of the box, place them on a flat and clean surface and identify each element of the kit (see content on the left).

SPECIMEN COLLECTION


Put the specimen collection tube (a) on the rack (marked position on the package).

Break off the extraction buffer tube (b).


Drop all the extraction buffer into the specimen collection tube (a).

Take the Specimen Bag (c) and open it. Cough several times and clear your throat deeply. Spit sputum (about 2 ml) into the Specimen Bag (c). Remove the dropper and use it to collect samples.

Place 5 drops (approx. 200ul) of the saliva into the sample collection tube. Cover the sample collection tube, and put the rack.

TEST PROCEDURE

Place the test device on a clean and horizontal surface. Reverse the specimen collection tube, extrude 3 drops of the prepared specimen into the specimen well (S) of the test cassette and start the timer.