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Hangzhou Immuno Biotech Co.,Ltd.

Hinzufügen: Room 401, 4F, Building 2, No. 28, No. 3 Street, Hangzhou Economic and Technological Development Zone, Hangzhou City, Zhejiang Province, China

Tel: 0571-85368996

E-mail: info@jgbiotech.com

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Lotus NL B.V.

Koningin Julianaplein 10,1e Verd, 2595AA, Den Haag, Niederlande.

E-mail: Peter@lotusnl.com

Importeur & Vertrieb:

Kagerlin GmbH

Basler Str. 7a

79540 Lörrach

Deutschland

www.kagerlin.com

Email: info@kagerlin.com


SARS-COV-2 Antigen Rapid Test COVID-19 Ag

SARS-CoV-2 Antigen Rapid Test (COVID-19) is a rapid chromatographic immunoassay for the qualitative detection of novel coronavirus SARS-CoV-2 in human nasal cavity.

Hangzhou Immuno Biotech Co.,Ltd.
Room401, 4F, Building 2, No.28, No.3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China

European Authorized Representative: Lotus NL B. V.
Koningin Julianaplein 10,1e Verd,2595AA,
The Hague, Netherlands.
E-mail: Peter@lotusnl.com

CONTENT OF THE KIT

PRECAUTIONS FOR USE


- For in vitro diagnostic use only. Do not use after expiration date.
- The Test strip should remain in the sealed pouch until use.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- The used Test Strip should be discarded according to national, state and local regulations.
- Humidity and temperature can adversely affect results.

INTERPRETATION OF RESULTS


- Positive (+): Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the T line region.
- *NOTE: The intensity of the colour in the test line regions may vary depending on the concentration of SARS-CoV-2 present in the specimen. Therefore, any shade of colour in the test line region should be considered positive and recorded as such.
- Negative (-): One coloured line appears in the control line region (C). No line appears in the T line region.
- Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

- PLEASE READ THE INSTRUCTIONS (and the information in the additional leaflet) THOROUGHLY BEFORE STARTING THE TEST PROCEDURE
- Have a watch or a timer (only equipment not provided).
- Take all the elements out of the box, place them on a flat and clean surface and identify each element of the kit (see content on the left).


SPECIMEN COLLECTION

TEST PROCEDURE

Bring the pouch (d) to room temperature before opening. Remove the test cassette (e) from the sealed pouch (d) and use it as soon as possible. Place the test cassette (e) on a clean and horizontal surface. Reverse the specimen collection tube (a), extrude 3 drops of the prepared specimen into the specimen well (S) of the test cassette (e) and start the timer. Wait for the colored line(s) to appear. Read results at 10 minutes. Do not interpret the result after 15 minutes.

IMMUNOBIO


SARS-COV-2 Antigen Rapid Test COVID-19 Ag

SARS-CoV-2 Antigen Rapid Test (COVID-19 Ag) is a rapid test for the qualitative detection of antigens to the novel coronavirus SARS-CoV-2 in sputum.

Hangzhou Immuno Biotech Co., Ltd.
Room 401, 4F, Building 2, No.28, No.3 Street, Hangzhou
Economic and Technological Development Zone, Hangzhou City,
Zhejiang Province, China

European Authorized Representative: Lotus NL B. V.
Koningin Julianaplein 10, 1e Verd, 2595AA,
The Hague, Netherlands.
E-mail: Peter@lotusnl.com

CONTENT OF THE KIT

PRECAUTIONS FOR USE

- For in vitro diagnostic use only. Do not use after expiration date.
- The Test strip should remain in the sealed pouch until use.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- The used Test Strip should be discarded according to national, state and local regulations.
- Humidity and temperature can adversely affect results.

INTERPRETATION OF RESULTS


- Positive (+): Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the T line region.
- *NOTE: The intensity of the colour in the test line regions may vary depending on the concentration of SARS-CoV-2 present in the specimen. Therefore, any shade of colour in the test line region should be considered positive and recorded as such.
- Negative (-): One coloured line appears in the control line region (C). No line appears in the T line region.
- Invalid: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

- PLEASE READ THE INSTRUCTIONS (and the information in the additional leaflet) THOROUGHLY BEFORE STARTING THE TEST PROCEDURE
- Have a watch or a timer (only equipment not provided).
- Take all the elements out of the box, place them on a flat and clean surface and identify each element of the kit (see content on the left).

SPECIMEN COLLECTION


Take the Specimen Bag (c) and open it. Cough several times and clear your throat deeply. Spit sputum (about 2 ml) into the Specimen Bag (c). Remove the dropper and use it to collect samples.

Place 5 drops (approx. 200ul) of the saliva into the sample collection tube. Cover the sample collection tube, and put the rack.

TEST PROCEDURE

Place the test device on a clean and horizontal surface. Reverse the specimen collection tube, extrude 3 drops of the prepared specimen into the specimen well (S) of the test cassette and start the timer.